
Product Datasheet

Mouse CD93 Protein (orb424647)

Description

Recombinant of mouse CD93 protein

Preservatives

CD93 protein solution (0.5mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.

Form/Appearance

Sterile filtered colorless solution.

Storage

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

Note

For research use only

Application notes

Recombinant & Natural Proteins

Protein Sequence

ADSQAVVCEG TACYTAHWGK LSAAEAQHRC
 NENGGNLATV KSEEEARHVQ QALTQLLTK
 APLEAKMGKF WIGLQREKGN CTYHDLPMRG
 FSWVGGGEDT AYSNWYKASK SSCIFKRCVS
 LILDLSLTPH PSHLPKWHEP PCGTPEAPGN
 SIEGFLCKFN FKGMCRLAL GGPGRVTTYT
 PFQATTSSLE AVPFASVANV ACGDEAKSET
 HYFLCNEKTP GIFHWGSSGP LCVSPKFGCS
 FNNGGCQQDC FEGGDGSFRC GCRPGFRLLD
 DLVTCASRNP CSSNPCTGGG MCHSVPLEN
 YTCRCPSGYQ LDSSQVHCVD IDECQDSPCA
 QDCVNTLGSF HCECWVGYQP SGPKEEACED
 VDECAAANSP CAQGCINTDG SFYCSCKEGY
 IVSGEDSTQC EDIDECS DAR GNPCDSLFCN
 TDGSFRCGCP PGWELAPNGV FCSRGTVFSE
 LPARPPQKED NDDRKESTMP PTEMPSSPSG
 SKDVSNRAQT TGLFVQSDIP TASVPLEIEI
 PSEVSDVWFE LGTYLPTTSG HSKPTHEDSV
 SAHSDDTDGQN LEHHHHHH

Purity

Greater than 95.0% as determined by SDS-PAGE.

Source

Sf9, Baculovirus cells

Expiration Date

12 months from date of receipt.