

Product Datasheet

Human CHMP2B Protein (orb424751)

Description	Recombinant of human CHMP2B protein
Preservatives	CHMP2B protein 1mg/ml is supplied in 20mM Tris-HCL, pH-8, 0.1M NaCl, 2mM DTT and 20% Glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SGLVPRGSH MASLFKKKTVD DDVIKEQNRE LRGTQRAIIR DRAALEKQEK QLELEIKKMA KIGNKEACKV LAKQLVHLRK QKTRTFVAVSS KVTSMSTQTK VMNSQMCMAG AMSTTAKTMQ AVNKKMDPQK TLQTMQNFQK ENMKMEMTEE MINDTLDDIF DGSDDEEESQ DIVNQLVDEI GIEISGKMAK APSAARSLPS ASTSKATISD EEIERQLKAL GVD
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.