

Product Datasheet

Human CLEC1B (Sf9) Protein (orb424790)

Description	Recombinant of human CLEC1B (Sf9) protein
Preservatives	CLEC1B protein solution (0.5mg/ml) containing Phosphate Buffered Saline (pH 7.4) and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	ADPSVMQRNY LQGENENRTG TLQQLAKRFC QYVVKQSELK GTFKGHKCSP CDTNWRYYGD SCYGFFRHNL TWEESKQYCT DMNATLLKID NRNIVEYIKA RTHLIRWVGL SRQKSNEVWK WEDGSVISEN MFEFLEDGKG NMNCAYFHNG KMHPTFCENK HYLMCERKAG MTKVDQLPHH HHHH
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells
Expiration Date	12 months from date of receipt.