

Product Datasheet

Human Cyclophilin E Protein (orb425000)

Description	Recombinant of human Cyclophilin E protein
Preservatives	Cyclophilin-E solution containing 20mM Tris pH-8.
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMTT KRVLYVGGLA EEVDDKVLHA AFIPFGDITD IQIPLDYETE KHRGFVFEV ELAEDAAAAI DNMNESELF RTIRVNLAKP MRIKEGSSRP VWSDDDWLKK FSGKLEENK EEEGSEPPKA ETQEGEPIAK KARSNPQVYM DIKIGNKPAG RIQMLLRSDV VPMTAENFRC LCTHEKGFGF KGSSFHRIIP QFMCQGGDFT NHNGTGGKSI YGKKFDDENF ILKHTGPGLL SMANSGPNTN GSQFFLTCDK TDWLDGKHVV FGEVTEGLDV LRQIEAQGSK DGKPKQKVII ADCGEYV
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.