

Product Datasheet

Human DHH (C23II) (His) Protein (orb425111)

Description	Recombinant of human DHH (C23II) (His) protein
Preservatives	The DHH (C23II) His solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 7.5), 0.15M NaCl, 1mM DTT and 10% glycerol.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMIIGPGR GPVGRRRYAR KQLVPLLYKQ FVPGVPERTL GASGPAEGRV ARGSERFRDL VPNYNPDIIF KDEENSGADR LMTERCKERV NALAIAVMNM WPGVRLRVTE GWDEDGHHAQ DSLHYEGRAL DITTSRDRN KYGLLARLAV EAGFDWVYYE SRNHVHVSVK ADNSLAVRAG G
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.