

Product Datasheet

Human DHH (C23II) Protein (orb425112)

Description	Recombinant of human DHH (C23II) protein
Preservatives	Lyophilized from a 0.2µm filtered concentrated solution in 1×PBS, pH 7.4.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Stability: Lyophilized DHH (C23II) although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution DHH (C23II) should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	IIGPGRGPVG RRRYARKQLV PLYKQFVPG VPERTLGASG PAEGRVARGS ERFRDLVPNY NPDIIFKDEE NSGADRLMTE RCKERVNALA IAVMNMWPGV RLRVTEGWDE DGHHAQDSLH YEGRALDITT SDRDRNKYGL LARLAVEAGF DWVYYESRNH VHVSVKADNS LAVRAGG
Purity	Greater than 95.0% as determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
Source	Escherichia Coli
Solubility (25°C)	It is recommended to reconstitute the lyophilized DHH (C23II) in sterile 18M-cm H ₂ O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	12 months from date of receipt.