

Product Datasheet

H3N2 Wyoming Protein (orb425956)

Description	Recombinant of H3N2 Wyoming protein
Preservatives	The Recombinant H3N2 A/Wyoming/2003/3 solution contains 10mM Sodium phosphate, pH 7.4, 150mM NaCl and 0.005% Tween-20.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: H3N2 A/Wyoming/2003/3 Recombinant should be stored at 4°C
Note	For research use only
Application notes	Viral Antigens- Influenza Virus
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	Baculovirus Insect Cells
Expiration Date	12 months from date of receipt.