

Product Datasheet

H7N7 Netherlands Protein (orb425959)

Description	Recombinant of H7N7 Netherlands protein
Preservatives	The Recombinant H7N7 A/Netherlands/219/03 solution contains 10mM Sodium phosphate, pH 7 and 150mM NaCl.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: H7N7 A/Netherlands/219/03 Recombinant should be stored at 4°C
Note	For research use only
Application notes	Immunological Activity: Western-Blot 0.1µg - 1µg per strip, ELISA 1µg/Well
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	Baculovirus Insect Cells
Expiration Date	12 months from date of receipt.