

Product Datasheet

HIV-1 gag p17-p24, gp41-gp120 Protein (orb426158)

Description	Recombinant of HIV-1 gag p17-p24, gp41-gp120 protein
Preservatives	20mM PBS pH 7.8, 20mM NaCl, 1mM DTT & 8M urea.
Form/Appearance	Sterile filtered colorless clear solution.
Storage	Stability: HIV-1 gag p17-p24, gp41-gp120 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles
Note	For research use only
Application notes	Applications: HIV-1 gag p17-p24, gp41-gp120 antigen is suitable for ELISA and Western blots, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems
Purity	Greater than 95.0% as determined by HPLC analysis and SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.