

Product Datasheet

Human IL 7 (His) Protein (orb426537)

Description	Recombinant of human IL 7 (His) protein
Preservatives	Interleukin -7 His is supplied in 1x PBS and 50% glycerol.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.