
Product Datasheet

Human IL23R Protein (orb426585)

Description

Recombinant of human IL23R protein

Preservatives

IL23R protein solution (0.25mg/ml) contains Phosphate Buffered Saline (pH7.4) and 10% glycerol.

Form/Appearance

Sterile filtered colorless solution.

Storage

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

Note

For research use only

Application notes

Cytokines And Growth Factors

Protein SequenceADPGITNINC SGHIWVEPAT IFKMGMNISI
YCQAAIKNCQ PRKLHFYKNG IKERFQITRI
NKTTARLWYK NFLEPHASMY CTAECPKHFO
ETLICGKDIS SGYPPDIPDE VTCVIYEYSG
NMTCTWNAGK LTYIDTKYVV HVKSLETEEE
QQYLTSSYIN ISTDSLQGGK KYLVWVQAAN
ALGMEESKQL QIHLDDIVIP SAAVISRAET
INATVPKTII YWDSQTTIEK VSCEMRYKAT
TNQTNVKEF DTNFTYVQQS EFYLEPNIKY
VFQVRCQETG KRYWQPWSSL FFHKTPETVP
QVTSKAFQHD TWNSGLTVAS ISTGHLTSDN
RGDIGLEPKS CDKTHTCPPC PAPELLGGPS
VFLFPPPKPD TLMISRTPEV TCVVVDVSHE
DPEVKFNWYV DGVEVHNAKT KPREEQYNST
YRVVSVLTVL HQDWLNGKEY KCKVSNKALP
APIEKTISKA KGQPREPQVY TLPPSRDELTA
KNQVSLTCLV KGFYPSDIAV EWESNGQPEN
NYKTTTPVLD SDGSFFLYSK LTVDKSRWQQ
GNVFSCSVMH EALHNHYTQK SLSLSPGKHH
HHHH**Purity**

Greater than 90.0% as determined by SDS-PAGE.

Source

Sf9, Baculovirus cells

Expiration Date

12 months from date of receipt.