

Product Datasheet

Human IL34 Protein (orb426593)

Description	Recombinant of human IL34 protein
Preservatives	The IL34 solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH8.0), 10% glycerol and 0.4M Urea.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMNE PLEMWPLTQN EECTVTGFLR DKLQYRSRLQ YMKHYFPINY KISVPYEGVF RIANVTRLQR AQSERELRY LWVLVLSLST ESVQDVLLEG HPSWKYLQEV QTLNINVOQG LTDVEVSPKV ESVLSLLNAP GPNLKLVRPK ALLDNCFRVM ELLYCSCCKQ SSVLNWQDCE VPSPQSCSPE PSLQYAATQL YPPPPWSPSS PPHSTGSRVP VRAQGEGLLP
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
Expiration Date	12 months from date of receipt.