

Product Datasheet

Human IRF 1 Protein (orb426662)

Description	Recombinant of human IRF 1 protein
Preservatives	1 mg/ml in 20mM Tris pH-8 and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Liquid IRF1 although stable at 10°C for 1 week, should be stored below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Please prevent freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MPITRMRMRP WLEMQINSNQ IPGLIWINK EMIFQIPWKHAAKHGWDINK DACLFRSWAI HTGRYKAGEK EPDPKTWKAN FRCAMNSLPD IEEVKDQSRN KGSSAVRVYR MLPP
Purity	Greater than 90.0% as determined by(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.