

Product Datasheet

Human LGALS7 (His) Protein (orb426863)

Description	Recombinant of human LGALS7 (His) protein
Preservatives	The Galectin-7 solution (1 mg/ml) 20mM Tris-HCl buffer (pH 8.0), 10% glycerol and 1mM DTT.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MSNVPHKSSL PEGIRPGTVL RIRGLVPPNA SRFHVNLLCG EEQGSDAALH FNPRLDTSEV VFNSKEQGSW GREERGPGVP FQRGQPFEVL IIASDDGFKA VVGDAQYHHF RHRLPLARVR LVEVGGDVQL DSVRIF
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.