

Product Datasheet

Human LZIC Protein (orb426950)

Description	Recombinant of human LZIC protein
Preservatives	The LZIC solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT and 10% glycerol.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SGLVPRGSH MGSMASRGKT ETSKLKQNL EQLDRLMQQL QDLEECREEL DTDEYEETKK ETLEQLSEFN DSLKKIMSGN MTLVDELSGM QLAIQAAISQ AFKTPEVIRL FAKKQPGQLR TRLAEMDRDL MVGKLERDLY TQQKVEILTA LRKLGEKLT DDEAFLSANA GAILSQFEKV STDLGSGDKI LALASFEVEK TTK
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.