

Product Datasheet

Human MAF1 Protein (orb426967)

Description	Recombinant of human MAF1 protein
Preservatives	The MAF1 solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0) and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHHH SSGLVPRGSH MKLLENSSE AINSQTLVET GDAHIIGRIE SYSCKMAGDD KHMFKQFCQE GQPHVLEALS PPQTSGLSPS RLSKSQGGEE EGPLSDKCSR KTLFYLIATL NESFRPDYDF STARSHEFSR EPSLSWVVNA VNCSLFSAVR EDFKDLKPQL WNAVDEEICL AECDIYSYNP DLSDPFGED GSLWSFNFF YNKRLKRIVF FSCRSISGST YTPSEAGNEL DMELGEEVE EESRSRGSGA EETSTMEEDR VPVICI
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
Expiration Date	12 months from date of receipt.