

Product Datasheet

Mouse Noggin Protein (orb427377)

Description	Mouse Noggin protein
Preservatives	Lyophilized from a 0.2µm filtered solution in 30% acetonitrile, 0.1% TFA.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Stability: Lyophilized Mouse Noggin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Mouse Noggin should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MQHYLHIRPAPSDNLPLVDLIEHPDPIFDPKEKDLNETLLRSLLGGHYDPGFMATSPPEDRP GGGGGPAGGAEDLAELDQLLRQRPSPGAMPSEIKGLEFSEGLAQGKKQRLSKLRRKLQM WLWSQTFPCPVLYAWNLDLGSRFWPRYVKVGSCFSKRSCSVPEGMVCKPSKSVHLTVLRWR CQRRGQRCGWIPYPIISECKCSC
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Solubility (25°C)	It is recommended to be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in 10mM HAc to a concentration of 0.1-1.0 mg/ml. Further dilutions should be made in appropriate buffered solutions.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)