

Product Datasheet

Human NOP16 Protein (orb427379)

Description	Recombinant of human NOP16 protein
Preservatives	NOP16 protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.2M NaCl, 50% glycerol, 2mM EDTA and 2mM DTT.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SGLVPRGSH MGSMPKAKGK TRRQKFGYSV NRKRLNRNAR RKAAPRIECS HIRHAWDHAK SVRQNLAEMG LAVDPNRAVP LRKRKVKAME VDIEERPKEK VRKPYVLNDL EAEASLPEKK GNTLSRDLID YVRYMVENHG EDYKAMARDE KNYQDTPKQ IRSKINVYKR FYPAEWQDFL DSLQKRKMEV E
Purity	Greater than 85.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.