

---

## Product Datasheet

### Mouse NT5E Protein (orb427419)

**Description**

Recombinant of mouse NT5E protein

**Preservatives**

NT5E protein solution (0.5mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.

**Form/Appearance**

Sterile Filtered colorless solution.

**Storage**

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

**Note**

For research use only

**Application notes**

Recombinant &amp; Natural Proteins

**Protein Sequence**

WELTILHTND VHSRLEQTS D DSTKCLNASL  
 CVGGVARLFT KVQIRKEEP NVLFLDAGDQ  
 YQGTIWFTVY KGLEVAHF MN ILGYDAMALG  
 NHEFDNGVEG LIDPLLRNVK FPILSANIKA  
 RGPLAHQISG LFLPSKVLSV GGEVVGIVGY  
 TSKETPFLSN PGTNLVFEDE ISALQPEVDK  
 LKTLNVNKII ALGHSGFEMD KLIAQKVRGV  
 DIVVGGHSNT FLYTGNPPSK EVPAGKY PFI  
 VTADDGRQVP VVQAYAFGKY LGYLKVEFDD  
 KGNVITSYGN PILLNSSIPE DATIKADINQ  
 WRIKLDNYST QELGRTIVYL DGSTQTCRFR  
 ECNMGNLICD AMINNNLRHP DEMFWNHVSM  
 CIVNGGGIRS PIDEKNGTI TWENLAAVLP  
 FGGTFDLVQL KGSTLKKAFE HSVHRYGQST  
 GEFLQVGGIH VVYDINRKPW NRVVQLEVLC  
 TKCRVPIYEP LEMDKVYKVT LPSYLANGGD  
 GFQMIKDELL KHDSGDQDIS VVSEYISKMK  
 VVYPAVEGRI KFSLEHHHHH H

**Purity**

Greater than 90% as determined by SDS-PAGE.

**Source**

Sf9, Baculovirus cells

**Expiration Date**

12 months from date of receipt.