

Product Datasheet

Human OTUB1 Protein (orb427503)

Description	Recombinant of human OTUB1 protein
Preservatives	The OTUB1 solution contains 20mM Tris buffer (pH 8.0) and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MAAEEPQQK QEPLGSDSEG VNCLAYDEAI MAQQDRIQQE IAVQNPLVSE RLELSVLYKE YAEDDNIYQQ KIKDLHKKYS YIRKTRPDGN CFYRAFGFSH LEALLDDSKE LQRFKAVSAK SKEDLVSQGF TEFTIEDFHN TFMDLIEQVE KQTSVADLLA SFNDQSTSDY LVVYLRLLLTS GYLQRESKFF EHFIEGGRTV KEFCQQEVEP MCKESDHIHI IALAQALSVS IQVEYMDRGE GGTTNPHIFP EGSEPKVYLL YRPGHYDILY K
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.