

Product Datasheet

Human PMP2 (His) Protein (orb427753)

Description	Recombinant of human PMP2 (His) protein
Preservatives	PMP2 His-Tag is supplied in 20mM Tris HCl pH-8 and 50% glycerol.
Form/Appearance	Sterile Filtered colorless liquid formulation.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.