

Product Datasheet

Mouse REN Protein (orb428120)

Description	Recombinant of mouse REN protein
Preservatives	REN protein solution (0.25mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	LPTRTATFER IPLKKMPSVR EILEERGVDM TRLSAEWGVF TKRPSLTNLT SPVVLTNLYN TQYYGEIGIG TPPQTFKVI DTGSANLWVP STKCSRLYLA CGIHSLYESS DSSSYMENG DFTIHYGSGR VKGFLSQDSV TVGGITVTQT FGEVTEPLI PFMLAKFDGV LGMGFPAQAV GGVTPVFDHI LSQGVLEEV FSVYYNRGSH LLGGEVVLGG SDPQHYQGNF HYVSISKTDS WQITMKGVSV GSSTLLCEEG CAVVVDTGSS FISAPTSSLK LIMQALGAKE KRIEEYVNC SQVPTLPDIS FDLGGRAYTL SSTDYVLQYP NRRDKLCTLA LHAMDIPPPT GPVWVLGATF IRKFYTEFDR HNNRIGFALA RHHHHHH
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells
Expiration Date	12 months from date of receipt.