

Product Datasheet

Human UBA5 Protein (orb429061)

Description	Recombinant of human UBA5 protein
Preservatives	The UBA5 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 10% glycerol and 50mM NaCl.
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SGLVPRGSH MGSMAESVE RLQQRVQELE RELAQERSLQ VPRSGDGGGG RVRIEKMSE VVDSNPYSRL MALKRMGIVS DYEKIRTFAY AIVGVGGVGS VTAEMLTRCG IGKLLFDYD KVELANMNRL FFQPHQAGLS KVQAAEHLR NINPDVLFV HNYNITTVENFQHFMDRISN GGLEEGKPV LVLSCVDNFE ARMTINTACN ELGQWME VSENAVSGHI QLIIPGESAC FACAPPLV ANIDEKTLKR EGVCAASLPT TMGVVAGIL QNVLKFLNFTVSYFLGYN AMQDFPTMS MKPNPQCDDR NCRKQEEYK KKVAALPKQE VIQEEEEIIH EDNEWGIELV SEVSEELKN FSGPVPDLPE GITVAYTIK QEDSVTELT VEDSGESLED LMAKMKNM
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.