

Product Datasheet

Human UBE2D4 Protein (orb429070)

Description	Recombinant of human UBE2D4 protein
Preservatives	The SRGN solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 2mM DTT, 0.1M NaCl and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MALKRIQKEL TDLQRDPPAQ CSAGPVGDDL FHWQATIMGP NDSPYQGGVF FLTIHFPTDY PFKPPKVAFT TKIYHPNINS NGSICLDILR SQWSPALTVS KVLISCSLL CDPNPDDPLV PEIAHTYKAD REKYNRLARE WTQKYAM
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
Expiration Date	12 months from date of receipt.