

Product Datasheet

Human UBE2J2 Protein (orb429079)

Description	Recombinant of human UBE2J2 protein
Preservatives	The UBE2J2 solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 200mM NaCl, 1mM DTT and 30% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMSSTSS KRAPTTATQR LKQDYLRICK DPVPYICAE LPSNILEWHY VVRGPEMTPY EGGYYHGKLI FPREFPFKPP SIYMITPNRGR FKCNTLRLCLS ITDFHPDTWN PAWSVSTILT GLLSFMVEKG PTLGSITSD FTKRQLAVQS LAFNLKDKVF CELFPEVVVE IKQKQKAQDE LSSRPQTLPL PDVVPDGETH LVQNGIQLLN GHAPGAVPNL AGLQQANRHH
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
Expiration Date	12 months from date of receipt.