

Product Datasheet

Human UBE2L6 Protein (orb429083)

Description	Recombinant of human UBE2L6 protein
Preservatives	The protein contains 20mM Tris pH-8, 0.1mM PMSF, 1mM DTT, and 10% glycerol.
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMASM RVVKELEDLQ KKPPPYLRNL SSDDANVLVW HALLLPDQPP YHLKAFNLRI SFPPEYFPKP PMIKFTTKIY HPNVDENGQI CLPIISENW KPCTKTCQVL EALNVLVNRN NIREPLRMDL ADLLTQNPPEL FRKNAEEFTL RFGVDRPS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.