

Product Datasheet

Human UBE2Q2 Protein (orb429086)

Description	Recombinant of human UBE2Q2 protein
Preservatives	UBE2Q2 protein solution (0.5mg/ml) containing Phosphate Buffer Saline (pH7.4), 30% glycerol and 1mM DTT.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SGLVPRGSH MGSMSVSGLK AELKFLASIF DKNHERFRIV SWKLDELHCQ FLVPQQGSPH SLPPPLTLHC NITESYPSSS PIWFDSEDP NLTSVLERLE DTKNNLLRQ QLKWLICELC SLYNLPKHL D VEMLDQPLPT GQNGTTEEVT SEEEEEEEM AEDIEDLDHY EMKEEEPISG KKSEDEGIEK ENLAILEKIR KTQRQDHLNG AVSGSVQASD RLMKELRDIY RSQSYKTGIY SVELINDSLY DWHVKLQKVD PDSPLHSDLQ ILKEKEGIEY ILLNFSFKDN FPDPPFVRV VLPVLSGGYV LGGGALCMEL LTKQGWSSAY SIESVIMQIN ATLVK GKARV QFGANKNQYN LARAQQSYNS IVQIHEKNGW YTPPKEDG
Purity	Greater than 85.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.