

Product Datasheet

Human UBE2R2 Protein (orb429087)

Description	Recombinant of human UBE2R2 protein
Preservatives	UBE2R2 protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 20% glycerol and 1mM DTT.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAQQQMT SSQKALMLEL KSLQEPEVEG FRITLVDESD LYNWEVAIFG PPNTLYEGGY FKAHIKFPID YPYSPPTFRF LTKMWHPNYI ENGDVCISIL HPPVDDPQSG ELPSEWNPNT QNVRTILLSV ISLLNEPNTF SPANVDASVM FRKWRDSKGG DKEYAEIIRK QVSATKAEAE KDGVKVPTTL AEYCIKTKVP SNDNSSLDDY DDLYDDDIDD EDEEEEDADC YDDDDSGNEE S
Purity	Greater than 85.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.