

## Product Datasheet

### Human UBE2T Protein (orb429089)

<b>Description</b>	Recombinant of human UBE2T protein
<b>Preservatives</b>	The UBE2T solution contains 20mM Tris-HCl pH-8, 1mM DTT, and 10% glycerol.
<b>Form/Appearance</b>	Sterile Filtered clear colorless solution.
<b>Storage</b>	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Enzymes
<b>Protein Sequence</b>	MQRASRLKRE LHMLATEPPP GITCWQDKDQ MDDLRAQILG GANTPYEKGV FKLEVIIPER YPFEPPIRF LTPIYHPNID SAGRICLDVL KLPPKGAWRP SLNIATVLTS IQLLMSEPNP DDPLMADISS EFKYNKPAFL KNARQWTEKH ARQKQKADEE EMLDNLPEAG DSRVHNSTQK RKASQLVGIE KKFHPDVLEH HHHHH
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	12 months from date of receipt.