

## Product Datasheet

### Human UBE2W Protein (orb429091)

<b>Description</b>	Recombinant of human UBE2W protein
<b>Preservatives</b>	UBE2W solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 10% glycerol, 2mM DTT and 0.1M NaCl.
<b>Form/Appearance</b>	Sterile filtered colorless solution.
<b>Storage</b>	Stability: UBE2W Human Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Enzymes
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MASMQRRLQK ELLALQNDPP PGMTLNEKSV QNSITQWIVD MEGAPGTLYE GEKFQLLFKF SSRYPFDSPQ VMFTGENIPV HPHVYSNGHI CLSILTEDWS PALSVQSVCL SIISMLSSCK EKRRPPDNSF YVRTCNKNPK KTKWWYHDDT C
<b>Purity</b>	Greater than 90.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	12 months from date of receipt.