

Product Datasheet

Human UFC1 Protein (orb429109)

Description	Recombinant of human UFC1 protein
Preservatives	UFC1 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: UFC1 Human Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MADEATTRVV SEIPVLKTNA GPRDRELWVQ RLKEEYQSLI RYVENNKNAD NDWFRLESNK EGTRWFGKCW YIHDLLKYEF DIEFDIPITY PTTAPEIAVP ELDGKTAKMY RGGKICLTDH FKPLWARNVP KFGLAHLMAL GLGPWLAVEI PDLIQGVIQ HKEKCNQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.