

Product Datasheet

Human UFD1L Protein (orb429110)

Description	Recombinant of human UFD1L protein
Preservatives	The UFD1L solution contains 20mM Tris-HCl buffer (pH 8.0), 0.1M NaCl, 1mM DTT and 30% glycerol.
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMFSFNMF DHPIPRVFQN RFSTQYRCFS VSMLAGPNDR SDVEKGGKII MPPSALDQLS RLNITYPMLF KLTNKNSDRM THCGVLEFVA DEGICYLPHW MMQNLLLEEG GLVQVESVNL QVATYSKFQP QSPDFLDITN PKAVLENALR NFACLTGTDV IAINYNEKIY ELRVMETKPD KAVSIECDM NVDFDAPLGY KEPERQVQHE ESTEGEADHS GYAGELGFRA FSGSGNRLDG KKKGVPESPS PIKPGDIKRG IPNYEFKLKG ITFIRNSRPL VKKVEEDEAG GRFVAFSGEG QSLRKKGRKP
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.