

Product Datasheet

Human USP46 Protein (orb429153)

| | |
|--------------------------|--|
| Description | Recombinant of human USP46 protein |
| Preservatives | The USP46 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 10% glycerol and 1mM DTT. |
| Form/Appearance | Sterile Filtered clear solution. |
| Storage | Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles |
| Note | For research use only |
| Application notes | Recombinant & Natural Proteins |
| Protein Sequence | MGSSHHHHHH SGLVPRGSH MGSMTVRNIA SICNMGTNAS ALEKDIGPEQ FPINEHYFGL VNFGNTCYCN SVLQALYFCR PFRENVLAYK AQKKKKNLL TCLADLFHSI ATQKKKVGVI PPKKFISRLR KENDLFDNYM QQDAHEFLNY LLNTIADILQ EEKKQEKQNG KLKNGNMNEP AENNKPELTW VHEIFQGTLT NETRCLNCET VSSKDEDFLD LSVDVEQNTS ITHCLRDFSN TETLCSEQKY YCETCCSKQE AQKRMVKKL PMILALHLKR FKYMEQLHRY TKLSYRVVFP LELRLFNTSS DAVNLDRMYD LVAVVVHCGS GPNRGHYITI VKSHGFWLLF DDDIVEKIDA QAIEEFYGLT SDISKNESEG YILFYQSRE |
| Purity | Greater than 90% as determined by SDS-PAGE. |
| Source | Escherichia Coli |
| Expiration Date | 12 months from date of receipt. |