

## Product Datasheet

### Human VPS29 Protein (orb429217)

<b>Description</b>	Recombinant of human VPS29 protein
<b>Preservatives</b>	The VPS29 solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 1mM DTT and 40% glycerol.
<b>Form/Appearance</b>	Sterile Filtered colorless solution.
<b>Storage</b>	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Recombinant & Natural Proteins
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSHMMLVLV LGDLHIPHRC NSLPAKFKKL LVPKIQHIL CTGNLCTKES YDYLKTLAGD VHIVRGDFDE NLNYPEQKVV TVGQFKIGLI HGHQVIPWGD MASLALLQRQ FDVDILISGH THKFEAFEHE NKFYINPGSA TGAYNALETN IIPSFVLMIDI QASTVVVTVYVY QLIGDDVKVE RIEYKKP
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Source</b>	E.coli
<b>Expiration Date</b>	12 months from date of receipt.