

## Product Datasheet

### Human VTI1B Protein (orb429227)

<b>Description</b>	Recombinant of human VTI1B protein
<b>Preservatives</b>	VTI1B protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 30% glycerol, 0.15M NaCl and 1mM DTT.
<b>Form/Appearance</b>	Sterile Filtered colorless solution.
<b>Storage</b>	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Recombinant & Natural Proteins
<b>Protein Sequence</b>	MGSSHHHHHH SGLVPRGSH MGSMASSAAS SEHFELKHEI FRGLHEDLQG VPERLLGTAG TEEKKLIRD FDEKQEQANE TLAEMEEELR YAPLSFRNPM MSKLRNYRKD LAKLHREVRS TPLTATPGGR GDMKYGIYAV ENEHMNRLQS QRAMLLQGTE SLNRATQSIE RSHRIATETD QIGSEIIEEL GEQRDQLERT KSRLVNTSEN LSKSRKILRS MSRKVTTNKL L
<b>Purity</b>	Greater than 85.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	12 months from date of receipt.