

Product Datasheet

Human YOD1 Protein (orb429269)

| | |
|--------------------------|---|
| Description | Recombinant of human YOD1 protein |
| Preservatives | The YOD1 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 30% glycerol and 1mM DTT. |
| Form/Appearance | Sterile Filtered colorless solution. |
| Storage | Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles |
| Note | For research use only |
| Application notes | Enzymes |
| Protein Sequence | MGSSHHHHHH SSGLVPRGSH MGSFMFGPAKG RHFGVHPAPG FPGGVSQQA A GTKAGPAGAW PVGSRDTMW RLRCKAKDGT HVLQGLSSRT RVRELQGQIA AITGIAPGGQ RILVGYPEEC LDLSNGDTIL EDLPIQSGDM LIIEEDQTRP RSSPAFTKRG ASSYVRETLV LTRTVVPAD NSCLFTSVYY VVEGGVLNPA CAPEMRRLIA QIVASDPDFY SEAILGKTNQ EYCDWIKRDD TWGGAIEISI LSKFYQCEIC VVDTQTVRID RFGEDAGYTK RVLLIYDGIH YDPLQRNFPD PDTPLTIFS SNDDIVLVQA LELADEARRR RQFTDVNRFT LRCMVCQKGL TGQAEAREHA KETGHTNFGV V |
| Purity | Greater than 90.0% as determined by SDS-PAGE. |
| Source | Escherichia Coli |
| Expiration Date | 12 months from date of receipt. |