

Product Datasheet

Human EBAG9 Protein (orb429368)

Description	Recombinant of human RCAS1 protein
Preservatives	The EBAG9 protein solution contains 20mM Tris-HCl pH-8, 1mM DTT, 2mM EDTA and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SGLVPRGSH MRSGRGRKLS GDQITLPTTV DYSSVPKQTD VEEWTSWDED APTSVKIEGG NGNVATQQNS LEQLEPDYFK DMTPTIRKTQ KIVIKKREPL NFGIPDGSTG FSSRLAATQD LPFIHQSSSEL GDLDTWQENT NAWEEEEEDAA WQAEVLRQQ KLADREKRAA EQQRKKMEKE AQRLMKKEQN KIGVKLS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.