
Product Datasheet

Human L1CAM Protein (orb429424)

Description

Recombinant of human CD106 protein

Preservatives

L1CAM protein solution (0.25mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.

Form/Appearance

Sterile filtered colorless solution.

Storage

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

Note

For research use only

Application notes

Recombinant & Natural Proteins

Protein Sequence

IQIPEELMEP PVITEQSPRR LVVFPDDIS
 LKCEASGKPE VQFRWTRDGV HFKPKEELGV
 TVYQSPHSGS FTITGNNSNF AQRFOGIYRC
 FASNKLGAM SHEIRLMAEG APKWPKETVK
 PVEVEEGESV VLPCNPPPSA EPLRIYWMNS
 KILHIKQDER VTMGQNGNLY FANVLTSDNH
 SDYICHAHFP GTRTIIQKEP IDLRVKATNS
 MIDRKPRLLF PTNSSSHLVA LQGQPLVLEC
 IAEGFPTPTI KWLRRPSGMP ADRVTYQNH
 KTLQLLKVGE EDDGEYRCLA ENSLGSARHA
 YYVTVEAAPY WLHKPQSHLY GPGETARLDC
 QVQGRPQPEV TWRINGIPVE ELAKDQKYRI
 QRGALILSNV QPSDTMVTQC EARNRHGLLL
 ANAYIYVVQL PAKILTADNQ TYMAVQGSTA
 YLLCKAFGAP VPSVQWLDED GTTVLQDERF
 FPYANGTLGI RDLQANDTGR YFCLAANDQN
 NVTIMANLKV KDATQITQGP RSTIEKKGSR
 VTFTCCASFD PSLQPSITWR GDGRDLQELG
 DSDKYFIEDG RLVHSLDYS DQGNYSVAS
 TELDVVESRA QLLVVGSPGP VPRLVLSDLH
 LLTQSQRVRS WSPAEDHNAP IEKYDIEFED
 KEMAPEKWYS LGKVPGNQTS TTLKLSPYVH
 YTFRVTAINK YGPGESPVS ETVVTPEAAP
 EKNPVDVKGE GNETTNMUIT WKPLRWMDWN
 APQVQYRVQW RPQGTGRPWQ EQIVSDPFLV
 VSNTSTFVPY EIKVQAVNSQ GKGPEPQVTI
 GYSGEDYPQA IPELEGIEIL NSSAVLVKWR
 PVDLAQVKGH LRGYNVTYWR EGSQRKHSKR
 HIIKDHVVVP ANTTSVILSG LRPYSSYHLE
 VQAFNGRGSG PASEFTFSTP EGVPGHPEAL
 HLECSNTSL LLRWQPPLSH NGVLTGYVLS
 YHPLDEGGKG QLSFNLRDPE LRTHNLTDLS
 PHLRYRFQLQ ATTKEGPGEA IVREGGTMAL
 SGISDFGNIS ATAGENYSVV SWVPKEGQCN
 FRFHILFKAL GEEKGGASLS POYVSYNOSS

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Source

Sf9, Baculovirus cells