

## Product Datasheet

### ADENOVIRUS antibody (Biotin) (orb429769)

<b>Description</b>	Goat polyclonal antibody against ADENOVIRUS conjugated to Biotin
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, ICC, IF, WB
<b>Immunogen</b>	Gradient purified, disrupted native adenovirus type 5 virions.
<b>Target</b>	ADENOVIRUS
<b>Preservatives</b>	Purified IgG conjugated to Biotin - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Biotin - liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.