

Product Datasheet

HEPATITIS C NS-5 ANTIGEN antibody (orb429782)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to HEPATITIS C NS-5 ANTIGEN |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IH-Fr, IHC-P, WB |
| Immunogen | Synthetic peptide corresponding to the NS-5 antigen region of the native molecule. |
| Target | HEPATITIS C NS-5 ANTIGEN |
| Preservatives | Purified IgG - liquid |
| Form/Appearance | Purified IgG - liquid |
| Concentration | 1.0 mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | Purified |
| Expiration Date | 12 months from date of receipt. |