

Product Datasheet

PARAINFLUENZA VIRUS 3 antibody (orb429851)

Description	Mouse monoclonal antibody to PARAINFLUENZA VIRUS 3
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Crude viral lysate, native.
Target	PARAINFLUENZA VIRUS 3
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.