
Product Datasheet

CGRP antibody (orb429900)

Description

Goat polyclonal antibody to CGRP

Species/Host

Goat

Reactivity

Rat

Conjugation

Unconjugated

Tested Applications

ELISA, ICC, IF, IHC-Fr, IHC-P

Immunogen

Synthetic rat Tyr-CGRP (23-37) conjugated to gamma globulin.

Target

CGRP

Preservatives

 Phosphate buffered saline. 0.09% Sodium Azide (NaN₃).

Form/Appearance

Purified Ig - liquid

Concentration

p>5.0 mg/ml

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity


Purified Ig prepared by affinity chromatography on Protein G

Uniprot ID
P01256
Dilution Range


ELISA: 1/500-2500

Expiration Date


12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Immunofluorescence analysis of mouse dor...



Immunofluorescence analysis of rat spina...