

## Product Datasheet

### MOUSE ANTI CORTISOL (orb429988)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CORTISOL                                                                                                   |
| <b>Species/Host</b>        | Mouse                                                                                                                                   |
| <b>Reactivity</b>          | Aves                                                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ELISA                                                                                                                                   |
| <b>Immunogen</b>           | Cortisol conjugated to bovine serum albumin.                                                                                            |
| <b>Target</b>              | CORTISOL                                                                                                                                |
| <b>Preservatives</b>       | Purified IgG - liquid                                                                                                                   |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                                                                                   |
| <b>Concentration</b>       | 1.0 mg/ml                                                                                                                               |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG1                                                                                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                                                                              |
| <b>Clone Number</b>        | BGN/C42                                                                                                                                 |
| <b>Purity</b>              | Purified                                                                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)