

Product Datasheet

CYTOMEGALOVIRUS LATE MATRIX ANTIGEN antibody (orb430000)

Description	Goat polyclonal antibody to CYTOMEGALOVIRUS LATE MATRIX ANTIGEN
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	Purified virions of strain AD169.
Target	CYTOMEGALOVIRUS LATE MATRIX ANTIGEN
Preservatives	Purified Ig - liquid
Form/Appearance	Purified Ig - liquid
Concentration	4 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P06473
Expiration Date	12 months from date of receipt.