

Product Datasheet

CYTOMEGALOVIRUS gB LATE CYTOPLASMIC ANTIGEN antibody (orb430003)

Description	Mouse monoclonal antibody to CYTOMEGALOVIRUS gB LATE CYTOPLASMIC ANTIGEN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Native cytomeglaovirus extracts from infected human cells..
Target	CYTOMEGALOVIRUS gB LATE CYTOPLASMIC ANTIGEN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.