

Product Datasheet

ENTEROVIRUS 71 antibody (orb430043)

Description	Mouse monoclonal antibody to ENTEROVIRUS 71
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FA, ICC, IF
Immunogen	Enterovirus 71 strain BrCr.
Target	ENTEROVIRUS 71
Preservatives	Ascitic Fluid - Liquid
Form/Appearance	Ascitic Fluid - Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Ascites
Expiration Date	12 months from date of receipt.