

Product Datasheet

FELINE IMMUNODEFICIENCY VIRUS p24 antibody (orb430069)

Description	Mouse monoclonal antibody to FELINE IMMUNODEFICIENCY VIRUS p24
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	Native, extracts
Target	FELINE IMMUNODEFICIENCY VIRUS p24
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Dilution Range	ELISA: 0.5 µg/well, FACS: 1:50/100, WB: 1:25/200
Expiration Date	12 months from date of receipt.