

Product Datasheet

HEPATITIS B CORE ANTIGEN antibody (orb430148)

Description	Mouse monoclonal antibody to HEPATITIS B CORE ANTIGEN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Purified hepatitis B core antigen.
Target	HEPATITIS B CORE ANTIGEN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.25 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.