

## Product Datasheet

# HEPATITIS B SURFACE ANTIGEN AD/AY antibody (orb430149)

<b>Description</b>	Goat polyclonal antibody to HEPATITIS B SURFACE ANTIGEN AD/AY
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Hepatitis B Surface Antigen Subtypes AD & AY
<b>Target</b>	HEPATITIS B SURFACE ANTIGEN AD/AY
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Concentration</b>	4.0 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Expiration Date</b>	12 months from date of receipt.