

Product Datasheet

HEPATITIS B RECEPTOR BINDING FRAGMENT PRE S antibody (orb430153)

Description	Rabbit polyclonal antibody to HEPATITIS B RECEPTOR BINDING FRAGMENT PRE S
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide corresponding to the amino acids 21-47 (PLGFFPDHQLDPAFGANSNNPDWDFNP(-C)) of the native molecule conjugated to Keyhole limpet Haemocyanin.
Target	HEPATITIS B RECEPTOR BINDING FRAGMENT PRE S
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.